STRÖME | Vitality Vibes Y3

Ohjekirja
Bruksanvisning
User manual
Brugsanvisning
Gebrauchsanleitung

Suomeksi

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
Safety

- Read the instructions and keep them for future reference. By following these safety instructions, you can use the device safely and prevent damage to the device or injury to persons.
- Before use, ensure that the package contains all parts listed in the technical specifications and that the device parts are intact. Do not use the device if any parts are missing or damaged.
- Do not attempt to repair the device parts yourself. Take the device to the retailer or an authorised service centre if any parts are faulty.
- Remove all protective plastics and stickers from the appliance.
- Use the appliance only in dry indoor environments at room temperature.
- Do not expose the appliance to water or other liquids.
- Do not expose the device to fire, direct sunlight, or other heat sources, excessive cold, or corrosive substances.
- Do not use this device instead of proper medical care. If you are in a medical treatment plan, talk to your doctor about the use of the device.
- Do not use the device to massage the head, intimate area or any body part recovering from surgery.
- The device will automatically shut down if the voltage is too low or the protection current too high, and after 15 minutes of continuous use.
- Charge the device fully whenever the battery level is low. Regularly charge the device even if you are not using it frequently.
- Charge the device in a dry area free from easily flammable or explosive substances, with a temperature between +5 °C and +40 °C.
- Ensure that the charging port, charger, and socket are dry before charging. Do not charge the device if the charging port, charger, or socket is wet.
- Charge the device's battery on a flat, stable surface. Do not let the device or charging cable hang over the edge of high surfaces during charging. Do not place the device or charging cable in walkways where they could pose a tripping hazard.
- Charge the device only with the supplied charger (model number Xiamen Xunheng Electronics Tech Co., Ltd XH0500-2000IUG).
- Plug the charger into a socket that is easily accessible.
- Do not leave the device charging unattended.
- If the device emits unusual smells or light during charging, stop charging immediately and contact service.
- Grip the plug when disconnecting the charging cable from the device or power source. Do not pull the charging cable.
- Disconnect the charging cable from the device and power source once the battery is fully charged.
- Never puncture the device's battery.
- The device can also be used by children over 8 years old and persons with physical or mental limitations, or persons without experience or knowledge of how to use the device, provided they have been instructed on its safe use or are under appropriate supervision and understand the risks associated with it. Children should not play with the appliance. Children should not clean or maintain the device without supervision.

Technical information

Model	Ströme Vitality Vibes Y3
Maximum power	60 W
Stroke length	6 mm
Speed	1800–3000 rpm
Battery	7.4 V Li-ion 1800 mAh / 13.32 Wh
Power consumption in off mode	0.2 W
Time after which the product automatically reaches off mode	0 min
Charging time	3,5 h
Operating time	approx. 4 h
Charging input	USB-C 5 V DC 2 A
Max noise level	50 dB
Dimensions	168.5 x 238.5 x 63 mm
Weight	680 g
Package contents	Massage gun, massage heads (6 pcs), charger, charging cable, user manual
Power supply	Xiamen Xunheng Electronics Tech Co., Ltd XH0500-2000IUG
Input voltage	100–240 V
Input frequency	50/60 Hz
Output voltage	5.0 V
Output current	2.0 A
Output power	10.0 W
Average active efficiency (230 V, 50 Hz)	79.77 %
No-load power consumption (230 V, 50 Hz)	0.06 W

Parts

1. Charging input
2. Indicator light (massage setting)
3. Indicator light (charging and power)
4. Power button
5. Massage head
6. Handle

Charging and use

Charging

Charge the device before first use and whenever the battery level is low. The battery level is low when the first power indicator light is on.

1. Plug the charging cable into the charging port and charger.
2. Plug the charger into a socket that is easily accessible. Charging begins and the power indicator light starts to blink.
3. All power indicator lights are on when the battery is fully charged. Disconnect the charging cable from the device and the charger. Disconnect the charger from the electric outlet as well.

How to use


Use the device to massage your back, arms, gluteal muscles, legs and feet.

- Insert the massage head directly into the holder. Do not change the massage head while the device is running.
- Hold the power button (4) to turn the device on and off.
- To change the massage setting, press the power button (4) quickly. The device indicates the massage setting in use with indicator lights. The device has six massage settings:


1. 1: Warms up muscles.
2. 2: Relaxes joints.
3. 3: Removes lactic acid.
4. 4: Deep massage.
5. 5: Professional massage.
6. 6: Automatic massage. The device detects the intensity of the massage and automatically changes the speed based on it.

- The device cannot be used while charging.
- The device shuts down automatically after 10 minutes of continuous use.
- The device produces sounds (e.g. vibration) during use. This is normal.


Tips


Massage back muscles three times for 10 to 20 seconds.


Massage arm muscles three times for 10 to 20 seconds.

Massage gluteal muscles four times for 20 to 60 seconds.

Massage leg muscles two times for 10 to 45 seconds.

Massage foot muscles three times for 20 to 30 seconds.

Maintenance, storage and recycling

- Clean the device with a soft dry cloth. Do not use other cleaning tools, detergents or sharp objects.
- Store the device only in dry indoor environments at room temperature.
- Store the device in its original package if you do not use it regularly.
- Safe disposal of the product: The appliance or its accessories must not be disposed of with regular household waste at the end of their life cycle. Take the appliance for recycling to the retailer or a local waste disposal facility responsible for recycling. For more information, contact your retailer or waste management company.

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